

Voice Recognition

Merlin Release Notes 2.3.52


What is it?

This feature uses the microphone connected to your device to recognise speech to transcribe into text.

This speech-to-text feature uses services available directly within your browser and is transmitted to services by the browser vendor. Please note that no audio spoken is transmitted to or stored within Merlin.


How do I use it?

When adding clinical notes, you will see a microphone button for each question.



The screenshot shows a window titled "Default Protocol" with two question sections. Each section has a "History" or "Examination" label with a microphone icon, a large text input area, and a "Select" dropdown menu with an "Answer" field and a "Clinical Code" field.

On first use, select 'Agree' to the disclaimer message.




The disclaimer dialog box contains a question mark icon and the following text: "This feature uses microphones connected to your device to recognise speech to transcribe into text. This speech-to-text feature uses services available directly within your browser and is transmitted to services by the browser vendor. No audio spoken is transmitted to or stored within Merlin. All transcribed text must be checked for accuracy; Vet Space Ltd take no responsibility or blame for the use of this feature. By selecting agree you confirm you understand this statement." Below the text are "Agree" and "Disagree" buttons.

Select the microphone icon to begin recording. Once selected, the microphone button will appear red to indicate voice recording is in process.



The screenshot shows the "Default Protocol" window with the microphone icon next to "Question 1" turned red, indicating that voice recording is in progress.

Once you have finished speaking, select the red microphone button to stop recording. Your speech will then be converted into text.



The screenshot shows a software window titled "Default Protocol". On the left, under "Question 1", there is a "History" section with a red microphone icon. Below it, a text input field contains the text "Patient has a skin condition". On the right, there is a "Select" section with an "Answer" dropdown menu and a "Clinical Code" input field.

Please check the text for accuracy before saving content into the Patient's history.