

PACS Integration

VetHub (PLH Medical)

Merlin 2.3.54





What is PACS Smart Link?

We have added an integration for VetHub working with PLH Medical; Smart Link allows you to access your PACS records where the Patient ID has been entered onto the Study records under the MRN field.

How does it work?

Accessing PACS via Smart Link

Users can view PACS from multiple ares of the system.

From the Reception Screen

Search for a Client, select a Patient and right-click on the Patient. From the menu, select 'View PACS'.



Treatment Menu

From the treatment screen, the 'View PACS' button will be avaiable.



Diary Screen and Day lists

From the Day lists or Diary screen, right click on an appointment. The 'View PACS' option will be available on this menu.





VetHub System

The view below is the view when logged into the VetHub system directly not via the Smart Link.

	Studies Activities 1 Analytics Administra		VetSpace VetSpace Integration-User * Sign C			
Studies VetSpace	Add Study • potie	ent name	Search advanced			
Patient Name Sec A ×	K MRN DOB DD-MM-YYYY	Accession #	Study Description	Referring Physician X Search Clear Save		
Study Stage Study Date All - All	Uploaded Date Starred T ~ All ~ All ~	hin Modality Tags All ~ All ~				
Actions -						
Patient	 Study 	Study Date	Uploaded	Actions		
SMALL ANIMAL 2 (M) DOB MRN 24855		16-09-2017 05:16 AM	10-07-2018 01:07 PM	💿 🔻 🖍 🔘 🗣 Reports 🕶 Unossigned 💌		

Filters can be used to locate your studies. Click the animal's name to view a study.

For example, 'Small Animal 2' can be viewed by clicking the eye icon or selecting the name.

For Smart Link to work you need to ensure that the MRN has field for the studies have the Patient ID from Merlin.

Selecting 'View PACS' will open a new tab and open the VetHub System and display any studies for the Patient.



If there are no studies for the patient selected, the web page will display:

Study not found.



How can I start using it?

If you are interested in this integration, please speak to your BDM or contact us on 0845 2000 725.

PACS is setup from the Administration menu. This will need to be configured for each site.

PACS Delete								
Cor	mpany	Туре	Url	Кеу	Site	Account Number	Account Name	
Cim	nar	AES	sdsd	•••••	Esher	343434343	dfgdfgdfgdfg	
Cim	nar	AES	tesst	•••••	Esher	11123344445	early dream	
Cim	nar	AES	https://vetspace.cimar.co.uk/api/v3/link/	•••••	Esher	a234727b-5261-4745-8121-9e1ef	CIMAR TEST	

The URL will be supplied from PLH Medical / VetHub (must contain /api/vs/link/ including trailing slash) e.g. <u>https://vetspace.cimar.co.uk/api/v3/link/</u> along with the API Key and Account Key. The account name is 'free text' and is used to describe the account you have added.

Once the above has been completed, choose the active PACS from the drop-down menu for Smart Link from Administration > Settings > PACS.

Categories	Site	Esher	~	H	Save	🕼 Reset			
Account	Dicture Archiving and Communication System								
Reporting	Picture Archiving and Communication System								
Diary	CI	MAR TEST			×				
🤱 Reception									
Treatment									
🧟 Referral									
🗱 Lab									
Vet Scan Lab									
Labelling									
SMS									
🕞 Email									
DayLists									
System									
Miscellaneous									
HTML Labellir									
Smart Flow									
Follow Up									
GDPR									
PACS									